Warum TESTA

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Why TESTA

TESTA Analytical Solutions is fully focused on Polymer Characterization and Particle Sizing Technologies, offering the most advanced and refined equipment for the most challenging research, development and QC Tasks. TESTA Analytical Solution's experts are known for their ability to help users solve the most intricate problems and applications, worldwide. Our know-how is the fundament of our success.

RELATED PRODUCTS

AYE dn/dc Refractometer

AYE Differential RI & Viscometer Combo

AYE Differential Refractive Index Detector

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*Deviations are possible and all items are subject to technical changes.

ReadyCal™ Kit*

Polymer science
at work

TESTA ANALYTICAL SOLUTIONS e.K.


TESTA
ANALYTICAL
SOLUTIONS

Certificate of Analysis

Product: ReadyCal-Kit PEO/PEG
 Part No: PSS-peokit1
 Lot No: peokit1-12

GPC/SEC - Calibration Curve


GPC/SEC - Calibration Table

Elution volume [ml]	Mp [Da]	Polymer Lot No:
20.84	969000	peo220414
21.83	517000	peo160114
23.08	222000	peo040808
24.95	86200	peo7032
26.33	42700	peo130810
27.71	18800	peg160-2
29.15	6690	peg60-3
30.73	1960	peg20-3
32.71	599	peg06-4
34.71	238	peg160113

Note:
 Mp = Molar mass at the peak maximum

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Why Standard Kit

GPC/SEC Chromatography is impossible without Polymer Standards. For QC Applications which mostly uses a single concentration detector, Polymer Standards are required for Calibration of the Columns which then allows calculation of Molecular Weights. More complex GPC/SEC Systems with Viscometers and MALS, require Polymer Standards for the determination of Instrument Constants and for System Validation. These are the facts which make GPC/SEC without Polymer Standards impossible.

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Why ReadyCal™ Kit*

A ReadyCal™ Kit* combines highest quality standards with easy and fast handling, reducing errors in sample preparation practically to zero. Three sets each of four different pre-mixed standards in 1.5 mL screw cap Vials build the strong backbone of a reliable calibration curve. Just add the desired amount of solvent, shake according to instruction and inject. Column calibration cannot be easier.